## Breast Shaper

A non-plastic surgery but effective breast enhancement  
 Gain a cup size in just 3 month

### Theory

Fats in the breast are floating due to gravity, they flow to our stomach at day time and to our back at night time while we sleep. The science is to use negative pressure to make space under the skin, sucking fats nearby to fill the space up. In addition, it keeps fats firm at its position, boosting circulation to the area, making it easier to grow fats as well as making the skin youthful.


### Applications

- 1 - For the purpose of enlargement, treatment should be made with a low negative pressure, and an 8 hour of treatment every day. Stop treatment for 1 or 2 days in every 7 days and continue again.
- 2 - A cup size could be gained with a period of 3 months. (e.g. from cup B to cup C)
- 3 - Try to gain some weight if you are very thin.
- 4 - For the purpose of firming, a 3 course of 2 hour treatment every week.
- 5 - To ensure breast development in youth, 2 hours of treatment every 2 days
- 6 - Apply lotion or baby oil on breast before and after treatment to reduce skin allergy (make sure you're not allergic to the lotion/baby oil you use)

### Introduction

- 1 - Lean forward slightly. Push surrounding fats towards breast while putting the first dome on. Each new dome comes with a skin adhesive on the edge to aid new application, the stickiness would disappear after about 10 times of use.
- 2 - Connect manual pump to the valve, press and release the pump twice.
- 3 - Repeats step 2 and 3 on the other side. Connect the 3 ways adaptor to both valves and the other end to the electric pump. Press bottom to have it operate automatically.
- 4 - Put the electric pump in the pocket inside the bra top. The electric pump functions by keeping appropriate pressure automatically. It operates while loosen and stops when the negative pressure is appropriate.
- 5 - Sufficient treatment would be made anytime for at least 2 hours with appropriate negative pressure.

### Cautions

- \* Over pressure at an extended period of time can cause skin problems.
- \* Do not use when getting skin diseases or allergy

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