## Medical beauty : Lift a little bit of wrinkles and scars everyday

It's natural that metabolism would slow down after puberty to reduce energy consumption, This can cause the slowing down of circulation. When the epidermis is not receiving enough essentials like oxygen and collagen due to poor circulation, the skin aim age at a mean quicker rate. Our Refine Pad are made of Elements of oxygen, carbon, hydrogen and silicone cross linked. It boost dermal circulation, encourages skin rejuvenate, self-hydrate and enhance skin elasticity. At the same time, it reduce excess water loss off, allowing enough nutrients to keep skin youthful.

- ✓ increase dermal circulation, retains moisture, boosts collagen and metabolism
- ✓ reduce and removes wrinkles, also prevents wrinkles from developing
- ✓ helps return skin to a youthful and supple appearance.
- ✓ reduce bag under eyes, avoids and reduces pimples
- ✓ Reduces effectively the appearance of scars
- ✓ 100% Medical instant muscle Silicone
- ✗ Additives or hormone

1. Clean skin thoroughly before applying, or apply over your general skin moisturizer could enhance its effect (oil-based substances may affect the adhesive quality)
2. Apply the pad to the desired area while sleeping or for more than 2 hours at any time. It is self-adhesive, viscosity may no longer exist after approximately 20 times of use. Patch drop off and stick on cloth or use on other skin care product could shorter it's life time.
3. Remove the pad slowly after use and put it back to the card inside the package.
4. it should be cleaned occasionally by gentle detergent and rinse with warm water so as to remove exfoliated skin cells, tissues and foreign objects.
5. Do not dry with a towel, may place back to card or use again without drying

1. Do not use when there is a skin problem or known is allergic to silicone.
2. Discontinue using it if a rash or irritation occurs.
3. Apply single use only to open wounds so as to prevent cross infections;
4. Make sure breathing is not interfered when using during sleeping.

Store in cool, ventilated and dry place, avoid contact with corrosives.  
Validity 1 years from date of manufacture.

EN ISO 13485:2012 Certificate No. Q1N 13 07 74603 006

SFDA Cert No. Guangdong 2011 2001

### Sign instructions:

Disposable

Sterile EO

Product number

Batch number

Expiry date

Date of manufacture

Name of manufacturer and Service :

EC Rep : Welkeng Ltd ([www.ce-marking.eu](http://www.ce-marking.eu))

29 Harley St., London W1G 9QR, UK

Mftr : Concord Medical Ltd

Address : Shop 9, Basement, Shing Lee Comm'l Bldg,  
No. 8 Wing Kut St., Sheung Wan, HK

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